

From: Separation Science <elearning.solutions@sepscience.com>
Sent: Thursday, December 22, 2011 1:17 PM
To: Hanchett, James (DPH)
Subject: Today in Separation Science - Latest Issue of HPLC Solutions

Dwell Differences

Online VersionForwardUnsubscribe

John Dolan's
Update

Separation Science e-Learning

HPLC Solutions

Issue 85
Dwell Differences

For the next few issues of *HPLC Solutions* we're going to look at the effect that the dwell volume has on gradient chromatograms. The dwell volume comprises all the HPLC system volume between the point the solvents are mixed to the head of the column. For high-pressure-mixing systems, this includes the mixer, connecting tubing, and sample loop. For low-pressure-mixing systems, we have the same components plus the volume of the pump heads and check valves. In the next *HPLC Solutions* we'll consider how to measure dwell volume, but first, let's look at some of the problems that can be created when the dwell volume differs between two HPLC system.
[Click for PDF>>](#)

Issue 84
Peak Purity

A reader recently asked a question about the use of the peak-purity function of his diode-array UV detector (DAD). The question related to whether or not he could detect the presence of enantiomers, stereoisomers, diastereoisomers, or epimers with the peak-purity function.
[Click for PDF>>](#)

Issue 83
Method Limits

A reader recently emailed me a question that went something like this: "I'm using LC-MS/MS to analyze a biomolecule that I have isolated from tissue by LC-MS/MS. The lowest point on the calibration curve that I prepared is 0.1 µM. All the samples we've tested have concentrations above 0.1 µM and the signal-to-noise (S/N) is 33 at this concentration. I know that S/N = 10 is considered the lower limit of quantification (LLOQ) – does this mean that I have to extend my curve downward, even if I'll never see concentrations that low?"
[Click for PDF>>](#)

Chromatography Forum' is now part of Separation Science. Click on the icons to read the latest 'hot topics'...

Liquid Chromatography

Gas Chromatography

Hyphenated Techniques

Sample Prep

Data Systems

ASK THE DOCTOR

If you have an analytical question for John Dolan...

CLICK HERE >>

FEATURED APPLICATIONS

DMABC Derivatization Agent Enables Improved Detection of Estradiol with ESI-MS
Company: Sigma-Aldrich
[Click to read>>](#)

Efficient SEC-UV-RI-MALS Analysis of Protein Aggregates
Company: Tosoh
[Click to read>>](#)

Rapid Analysis of Polycyclic Aromatic Hydrocarbons
Company: Knauer
[Click to read>>](#)

Published by Eclipse Business Media Ltd
Frederick House | Princes Court | Beam Heath Way | Nantwich | Cheshire CW5 6PQ | United Kingdom
20 Maxwell Road | #09-17 Maxwell House | Singapore 069113
Copyright © 2011 Eclipse Business Media Ltd. All rights reserved.

To sponsor or advertise in this e-newsletter contact dean.grimes@sepscience.com
For editorial contributions, contact david.hills@sepscience.com

This message was sent from Separation Science to james.hanchett@state.ma.us. It was sent from: Eclipse Business Media Ltd, Frederick House, Princes Court, Beam Heath Way, Nantwich, Cheshire CW5 6PQ, United Kingdom. You can modify/update your subscription via the link below.

[Unsubscribe](#) | [Forward To a Friend](#)

1

Folk_OIG_PRR_093999